

TACR3 antibody



Catalog Number: 115972

Product name

TACR3 antibody

Specificity

Human; other species not tested.

Antibody description

TACR3 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue.

Preparation

This antibody was obtained by immunization of TACR3 recombinant protein (Accession Number: NM_001059). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:115972(TACR3 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:115972(TACR3 Antibody) at dilution of 1:50