

# TACR3 antibody

Catalog Number: 115972

---

## Product name

TACR3 antibody

## Specificity

Human; other species not tested.

## Antibody description

TACR3 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue.

## Preparation

This antibody was obtained by immunization of TACR3 recombinant protein (Accession Number: NM\_001059). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:115972(TACR3 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:115972(TACR3 Antibody) at dilution of 1:50