

TADA1L antibody

Catalog Number: 115973

Product name

TADA1L antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

TADA1L Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, NIH/3T3 cells. Positive IHC detected in human testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 38 kDa




Preparation

This antibody was obtained by immunization of TADA1L recombinant protein (Accession Number: NM_053053). Purification method: Antigen affinity purified.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115973(TADA1L Antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115973(TADA1L Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

Immunohistochemistry of paraffin-embedded

TADA1L antibody

Catalog Number: 115973

human testis tissue slide using Catalog No:115973(TADA1L Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:115973(TADA1L Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)