

TADA3L antibody



Catalog Number: 115976

Product name

TADA3L antibody

Immunogen

[Human TADA3L Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TADA3L Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human brain tissue, NIH/3T3 cells. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of TADA3L recombinant protein (Accession Number: NM_133481). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115976(TADA3L antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:115976(TADA3L antibody) at dilution of 1:100 (under 10x lens)