

# TAF12 antibody

Catalog Number: 115978

## Product name

TAF12 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TAF12 Rabbit Polyclonal antibody. Positive IP detected in Y79 cells. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of TAF12 recombinant protein (Accession Number: NM\_001135218). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:115978(TAF12 antibody) at dilution of 1:300



IP Result of anti-TAF12 (IP:Catalog No:115978, 4ug; Detection:Catalog No:115978 1:300) with Y79 cells lysate 2000ug.