

# TAF1C antibody



Catalog Number: 115983

## Product name

TAF1C antibody

## Specificity

human; other species not tested.

## Antibody description

TAF1C Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells. Observed molecular weight by Western-blot: 95kd

## Preparation

This antibody was obtained by immunization of TAF1C recombinant protein (Accession Number: NM\_139353). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:115983(TAF1C antibody) at dilution of 1:300