

TAF1L antibody

Catalog Number: 115984

Product name

TAF1L antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TAF1L Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, DU 145 cells. Observed molecular weight by Western-blot: 250-260 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_153809). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115984(TAF1L antibody) at dilution of 1:800