

# TAF3 antibody



Catalog Number: 115986

## Product name

TAF3 antibody

## Specificity

Human; other species not tested.

## Antibody description

TAF3 Rabbit Polyclonal antibody. Positive WB detected in cells. Observed molecular weight by Western-blot: 140 kDa

## Preparation

This antibody was obtained by immunization of TAF3 recombinant protein (Accession Number: BC024661). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: N/A

## Validations



WB results of TAF3 (Catalog No:115986) with siRNA sample. (140 kDa band absent in siRNA sample)