

TAF5L antibody



Catalog Number: 115987

Product name

TAF5L antibody

Immunogen

[Human TAF5L Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

TAF5L Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of TAF5L recombinant protein (Accession Number: NM_001025247). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115987(TAF5L antibody) at dilution of 1:500