

TAF9B antibody

Catalog Number: 115991

Product name

TAF9B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TAF9B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive IHC detected in human breast cancer tissue.

Preparation

This antibody was obtained by immunization of TAF9B recombinant protein (Accession Number: NM_015975). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP


Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115991(TAF9B antibody) at dilution of 1:100



(under 10x lens)

IP Result of anti-TAF9B (IP:Catalog No:115991, 4ug; Detection:Catalog No:115991 1:300) with HeLa cells lysate 4000ug.