

TFAP2E antibody



Catalog Number: 116002

Product name

TFAP2E antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TFAP2E Rabbit Polyclonal antibody. Positive WB detected in A375 cells, mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of TFAP2E recombinant protein (Accession Number: BC041175). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:116002(TFAP2E Antibody) at dilution of 1:600