

TFAP4 antibody



Catalog Number: 116003

Product name

TFAP4 antibody

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Specificity

human, mouse, rat; other species not tested.

Validations

Antibody description

TFAP4 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 39kd




Preparation

This antibody was obtained by immunization of TFAP4 recombinant protein (Accession Number: NM_003223). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116003(TFAP4 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116003(TFAP4 antibody) at dilution of 1:50 (under 10x lens)

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications