

# TFAP4 antibody



Catalog Number: 116003

## Product name

TFAP4 antibody

ELISA, WB, IHC

## Specificity

human, mouse, rat; other species not tested.

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


## Antibody description

TFAP4 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 39kd

## Validations

## Preparation

This antibody was obtained by immunization of TFAP4 recombinant protein (Accession Number: NM\_003223). Purification method: Antigen affinity purified.



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116003(TFAP4 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116003(TFAP4 antibody) at dilution of 1:50 (under 10x lens)