

TFB2M antibody



Catalog Number: 116004

Product name

TFB2M antibody

Immunogen

[Human TFB2M Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TFB2M Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse heart tissue. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of TFB2M recombinant protein (Accession Number: NM_022366). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116004(TFB2M Antibody) at dilution of 1:600