

# TFDP1 antibody

Catalog Number: 116006

---

## Product name

TFDP1 antibody

## Immunogen

[Human TFDP1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TFDP1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue.

## Preparation

This antibody was obtained by immunization of TFDP1 recombinant protein (Accession Number: XM\_047430563). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116006(TFDP1 antibody) at dilution of 1:50 (under 10x lens)