

TFF2 antibody

Catalog Number: 116011

Product name

TFF2 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TFF2 Rabbit Polyclonal antibody. Positive IHC detected in human stomach cancer tissue, human stomach tissue, mouse liver tissue, mouse pancreas tissue. Positive WB detected in human stomach tissue. Observed molecular weight by Western-blot: 18-20 kDa

Validations




Preparation

This antibody was obtained by immunization of TFF2 recombinant protein (Accession Number: NM_005423). Purification method: Antigen affinity purified.

human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:116011(TFF2 antibody) at dilution of 1:100

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Immunohistochemical of paraffin-embedded human stomach cancer using Catalog No:116011(TFF2 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human stomach cancer using Catalog

TFF2 antibody



Catalog Number: 116011

No:116011(TFF2 antibody) at dilution of 1:50 (under 10x lens)