

TFPI antibody



Catalog Number: 116015

Product name

TFPI antibody

Specificity

Human; other species not tested.

Antibody description

TFPI Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of TFPI recombinant protein (Accession Number: NM_001032281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116015(TFPI antibody) at dilution of 1:500



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:116015(TFPI Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)