

TGF-alpha antibody



Catalog Number: 116018

Product name

TGF-alpha antibody

Specificity

Human; other species not tested.

Antibody description

TGF-alpha Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive IHC detected in human colon tissue, human breast cancer tissue, human skin tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003236). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemistry of paraffin-embedded human colon slide using Catalog No:116018(TGFA Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human colon slide using Catalog No:116018(TGFA Antibody) at dilution of 1:50



Immunofluorescent analysis of (-20oc Ethanol) fixed A431 cells using Catalog No:116018(TGFA Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)