

# TGFBR3-Specific antibody



Catalog Number: 116024

## Product name

TGFBR3-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

TGFBR3-Specific Rabbit Polyclonal antibody. Positive FC detected in A549 cells. Positive WB detected in A549 cells, A431 cells. Observed molecular weight by Western-blot: 110-150 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003243). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116024(TGFBR3-Specific antibody) at dilution of 1:300



1X10<sup>6</sup> A549 cells were stained with 0.2ug TGFBR3-Specific antibody (Catalog No:116024, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.