

TGIF2 antibody



Catalog Number: 116026

Product name

TGIF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TGIF2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of TGIF2 recombinant protein (Accession Number: NM_001199515). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



Immunofluorescent analysis of HepG2 cells, using TGIF2 antibody Catalog No:116026 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116026(TGIF2 Antibody) at dilution of 1:600