

TH1L antibody



Catalog Number: 116035

Product name

TH1L antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

TH1L Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 66 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of TH1L recombinant protein (Accession Number: NM_016397). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116035(TH1L antibody) at dilution of 1:500

Clonality