

THRA antibody



Catalog Number: 116052

Product name

THRA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

THRA Rabbit Polyclonal antibody. Positive WB detected in A431 cells, mouse brain tissue. Observed molecular weight by Western-blot: 58kd, 50kd

Preparation

This antibody was obtained by immunization of THRA recombinant protein (Accession Number: NM_003250). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:116052(THRA antibody) at dilution of 1:2000