

# Thrombospondin 3-Specific antibody



Catalog Number: 116055

## Product name

Thrombospondin 3-Specific antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Thrombospondin 3-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 40-45 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_007112).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116055(THBS3-Specific antibody) at dilution



of 1:1000

IP Result of anti-THBS3-Specific (IP:Catalog No:116055, 4ug; Detection:Catalog No:116055 1:1000) with mouse kidney tissue lysate 4000ug.