

THRSP antibody

Catalog Number: 116057

Product name

THRSP antibody

Immunogen

[Human THRSP Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

THRSP Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, Transfected HEK-293 cells. Positive IF detected in MCF-7 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 20kd

Preparation

This antibody was obtained by immunization of THRSP recombinant protein (Accession Number: NM_003251). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116057(THRSP Antibody) at dilution of 1:1000



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116057(THRSP Antibody) at dilution of 1:700




THRSP antibody



Catalog Number: 116057

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116057(THRSP Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells using Catalog No:116057(THRSP Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)