

TM9SF4 antibody



Catalog Number: 116080

Product name

TM9SF4 antibody

Specificity

Human; other species not tested.

Antibody description

TM9SF4 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, A375 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of TM9SF4 recombinant protein (Accession Number: NM_001363731). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:300-1:1200

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:116080(TM9SF4 Antibody) at dilution of 1:1000