

# TM9SF4 antibody



Catalog Number: 116080

## Product name

TM9SF4 antibody

## Specificity

Human; other species not tested.

## Antibody description

TM9SF4 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, A375 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of TM9SF4 recombinant protein (Accession Number: NM\_001363731). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:300-1:1200

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:116080(TM9SF4 Antibody) at dilution of 1:1000