

TMCC2 antibody



Catalog Number: 116083

Product name

TMCC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMCC2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 85-90 kDa

Preparation

This antibody was obtained by immunization of TMCC2 recombinant protein (Accession Number: NM_014858). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116083(TMCC2 Antibody) at dilution of 1:1000



-	+	-	+	Overexpression of TMCC2-his-myc
+	+	-	-	Rabbit Anti human TMCC2 polyclonal antibody
-	-	+	+	Mouse Anti his-tag monoclonal antibody

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116083(TMCC2 Antibody) at dilution of 1:1000