

# TMEM11 antibody



Catalog Number: 116093

## Product name

TMEM11 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

TMEM11 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse placenta tissue. Observed molecular weight by Western-blot: 22 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of TMEM11 recombinant protein (Accession Number: NM\_003876). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116093(TMEM11 antibody) at dilution of 1:1000

## Clonality