

# TMEM176A antibody



Catalog Number: 116107

## Product name

TMEM176A antibody

## Specificity

Human; other species not tested.

## Antibody description

TMEM176A Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of TMEM176A recombinant protein (Accession Number: NM\_018487). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of Hela cells, using TMEM176A antibody Catalog No:116107 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116107(TMEM176A Antibody) at dilution of 1:700