

TIMM44 antibody



Catalog Number: 116133

Product name

TIMM44 antibody

Immunogen

[Human TIMM44 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TIMM44 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, Jurkat cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of TIMM44 recombinant protein (Accession Number: NM_006351). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116133(TIMM44 antibody) at dilution of 1:500



Immunofluorescent analysis of A549 cells, using TIMM44 antibody Catalog No:116133 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).