

# TIMM8A antibody

Catalog Number: 116135

## Product name

TIMM8A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TIMM8A Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human brain tissue, human liver tissue, mouse testis tissue, rat kidney tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 11 kDa

## Preparation

This antibody was obtained by immunization of TIMM8A recombinant protein (Accession Number: NM\_004085). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116135(TIMM8A antibody) at dilution of 1:500



IP Result of anti-TIMM8A (IP:Catalog No:116135, 4ug; Detection:Catalog No:116135 1:500) with mouse testis tissue lysate 4000ug.