

TIMP1 antibody



Catalog Number: 116138

Product name

TIMP1 antibody

Specificity

Human; other species not tested.

Antibody description

TIMP1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human prostate cancer tissue. Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of TIMP1 recombinant protein (Accession Number: NM_003254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200


Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116138(TIMP1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116138(TIMP1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:116138(TIMP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure

TIMP1 antibody



Catalog Number: 116138

Goat Anti-Rabbit IgG(H+L)