

TIMP1 antibody

Catalog Number: 116139

Product name

TIMP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TIMP1 Rabbit Polyclonal antibody. Positive IF detected in BxPC-3 cells. Positive IHC detected in human stomach cancer tissue, human pancreas tissue, human thyroid cancer tissue.

Preparation

This antibody was obtained by immunization of TIMP1 recombinant protein (Accession Number: BC007097). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using Catalog No:116139(TIMP1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using Catalog No:116139(TIMP1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of BxPC-3 cells using Catalog No:116139(TIMP1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)