

Timp-3 antibody

Catalog Number: 116141

Product name

Timp-3 antibody

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Timp-3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human kidney tissue, human lung cancer tissue, human oesophagus cancer tissue, human stomach tissue.


Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116141 (Timp-3 Antibody) at dilution of 1:50 (under 10x lens)

Preparation

This antibody was obtained by immunization of Timp-3 recombinant protein (Accession Number: NM_000362). Purification method: Antigen affinity purified.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116141 (Timp-3 Antibody) at dilution of 1:50 (under 40x lens)

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications