

TNNT1 antibody



Catalog Number: 116154

Product name

TNNT1 antibody

Polyclonal

Ig Type

Specificity

Human, Mouse, Rat; other species not tested.

Rabbit IgG

Applications

Antibody description

TNNT1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 30-33kd

ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of TNNT1 recombinant protein (Accession Number: BC022086). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116154(TNNT1 antibody) at dilution of 1:1500

Clonality