

TNS1 antibody



Catalog Number: 116159

Product name

TNS1 antibody

Specificity

Human; other species not tested.

Antibody description

TNS1 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 186 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_022648). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116159(TNS1 Antibody) at dilution of 1:300