

TMEM201 antibody



Catalog Number: 116175

Product name

TMEM201 antibody

Immunogen

[Human TMEM201 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TMEM201 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 43-45 kDa

Preparation

This antibody was obtained by immunization of TMEM201 recombinant protein (Accession Number: XM_011540910). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116175(TMEM201 Antibody) at dilution of 1:100