

# TMEM24 antibody



Catalog Number: 116180

## Product name

TMEM24 antibody

## Immunogen

[Human C2CD2L Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

TMEM24 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 80-100 kDa

## Preparation

This antibody was obtained by immunization of TMEM24 recombinant protein (Accession Number: NM\_001382613). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116180(C2CD2L Antibody) at dilution of 1:1000