

TMEM41A antibody



Catalog Number: 116188

Product name

TMEM41A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TMEM41A Rabbit Polyclonal antibody. Positive WB detected in mouse small intestine tissue. Observed molecular weight by Western-blot: 26-28kd

Preparation

This antibody was obtained by immunization of TMEM41A recombinant protein (Accession Number: NM_080652). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:116188(TMEM41A antibody) at dilution of 1:800