

TMEM55B antibody



Catalog Number: 116192

Product name

TMEM55B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TMEM55B Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of TMEM55B recombinant protein (Accession Number: NM_144568). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116192(TMEM55B Antibody) at dilution of 1:300