

TPSAB1 antibody



Catalog Number: 116217

Product name

TPSAB1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TPSAB1 Rabbit Polyclonal antibody. Positive IHC detected in human appendix tissue, human tonsillitis tissue. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse colon tissue. Observed molecular weight by Western-blot: 32-38 kDa

Validations




Preparation

This antibody was obtained by immunization of TPSAB1 recombinant protein (Accession Number: BC028059). Purification method: Antigen affinity purified.

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116217(TPSAB1 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human appendix using Catalog No:116217(TPSAB1 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemical of paraffin-embedded

TPSAB1 antibody

Catalog Number: 116217

human appendix using Catalog
No:116217(TPSAB1 antibody) at dilution of 1:100

IP Result of anti-TPSAB1 (IP:Catalog No:116217,
4ug; Detection:Catalog No:116217 1:800) with
mouse liver tissue lysate 4000ug.

