

TPSD1 antibody

Catalog Number: 116218

Product name

TPSD1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TPSD1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive IP detected in A549 cells. Observed molecular weight by Western-blot: 26-30 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_012217). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116218(TPSD1 Antibody) at dilution of 1:300



IP Result of anti-TPSD1 (IP:Catalog No:116218, 3ug; Detection:Catalog No:116218 1:300) with A549 cells lysate 920ug.