

TRADD antibody



Catalog Number: 116228

Product name

TRADD antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

TRADD Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, mouse kidney tissue. Observed molecular weight by Western-blot: 34 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of TRADD recombinant protein (Accession Number: NM_003789). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116228(TRADD antibody) at dilution of 1:500

Clonality