

# TRAF1 antibody



Catalog Number: 116229

## Product name

Polyclonal

TRAF1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

TRAF1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 46-50 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of TRAF1 recombinant protein (Accession Number: NM\_005658). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116229(TRAF1 antibody) at dilution of 1:500