

TOMM7 antibody

Catalog Number: 116243

Product name

TOMM7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TOMM7 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human gliomas tissue, human brain tissue. Observed molecular weight by Western-blot: 6 kDa

Preparation

This antibody was obtained by immunization of TOMM7 recombinant protein (Accession Number: NM_019059). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116243(TOMM7 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116243(TOMM7 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116243(TOMM7 antibody) at dilution of 1:50 (under 40x lens)