

TOR1AIP2 antibody



Catalog Number: 116252

Product name

TOR1AIP2 antibody

Specificity

Human; other species not tested.

Antibody description

TOR1AIP2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of TOR1AIP2 recombinant protein (Accession Number: NM_001349934). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:

WB: 1:200-1:1000


IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116252(TOR1AIP2 Antibody) at dilution of

1:300



Immunofluorescent analysis of A431 cells using Catalog No:116252(TOR1AIP2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)