

# TORC1 antibody



Catalog Number: 116255

## Product name

TORC1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TORC1 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue, A431 cells, Jurkat cells, L02 cells, mouse brain tissue. Observed molecular weight by Western-blot: 75 kDa

## Preparation

This antibody was obtained by immunization of TORC1 recombinant protein (Accession Number: NM\_015321). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116255(TORC1 antibody) at dilution of 1:800



IP Result of anti-TORC1 (IP:Catalog No:116255, 4ug; Detection:Catalog No:116255 1:800) with mouse kidney tissue lysate 4000ug.