

TOX2 antibody

Catalog Number: 116257

Product name

TOX2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TOX2 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of TOX2 recombinant protein (Accession Number: BC007636). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116257(TOX2 Antibody) at dilution of 1:300



IP Result of anti-TOX2 (IP:Catalog No:116257, 4ug; Detection:Catalog No:116257 1:300) with mouse liver tissue lysate 4000ug.