

TPCN1 antibody



Catalog Number: 116261

Product name

TPCN1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TPCN1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 94 kDa

Preparation

This antibody was obtained by immunization of TPCN1 recombinant protein (Accession Number: NM_017901). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116261(TPCN1 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:116261(TPCN1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)