

# TMOD1 antibody



Catalog Number: 116269

## Product name

TMOD1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TMOD1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human brain tissue, mouse brain tissue, mouse heart tissue. Observed molecular weight by Western-blot: 41kd

## Preparation

This antibody was obtained by immunization of TMOD1 recombinant protein (Accession Number: NM\_001166116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116269(TM0D1 Antibody) at dilution of 1:600